

Catalogue numbers for GORE® EXCLUDER® AAA Endoprosthesis

Trunk - Ipsilateral Leg Endoprosthesis

GORE® C3® Delivery System Catalogue Number	Aortic Endoprosthesis Diameter (mm)	Iliac Endoprosthesis Diameter (mm)	Endoprosthesis Length (cm)
RLT231212	23	12	12
RLT231214	23	12	14
RLT231216	23	12	16
RLT231218	23	12	18
RLT231412	23	14.5	12
RLT231414	23	14.5	14
RLT231416	23	14.5	16
RLT231418	23	14.5	18
RLT261212	26	12	12
RLT261214	26	12	14
RLT261216	26	12	16
RLT261218	26	12	18
RLT261412	26	14.5	12
RLT261414	26	14.5	14
RLT261416	26	14.5	16
RLT261418	26	14.5	18
RLT281212	28.5	12	12
RLT281214	28.5	12	14
RLT281216	28.5	12	16
RLT281218	28.5	12	18
RLT281412	28.5	14.5	12
RLT281414	28.5	14.5	14

GORE® C3® Delivery System Catalogue Number	Aortic Endoprosthesis Diameter (mm)	Iliac Endoprosthesis Diameter (mm)	Endoprosthesis Length (cm)
RLT281416	28.5	14.5	16
RLT281418	28.5	14.5	18
RLT311413	31	14.5	13
RLT311415	31	14.5	15
RLT311417	31	14.5	17
RLT351414	35	14.5	14
RLT351416	35	14.5	16
RLT351418	35	14.5	18

Contralateral Leg Endoprosthesis*

Catalogue Number	Distal Diameter (mm)	Endoprosthesis Length (cm)
PLC121000	12	10
PLC121200	12	12
PLC121400	12	14
PLC141000	14.5	10
PLC141200	14.5	12
PLC141400	14.5	14
PLC161000*	16	9.5
PLC161200*	16	11.5
PLC161400*	16	13.5
PLC181000*	18	9.5
PLC181200*	18	11.5
PLC181400*	18	13.5
PLC201000*	20	9.5
PLC201200*	20	11.5
PLC201400*	20	13.5
PLC231000*	23	10

Catalogue Number	Distal Diameter (mm)	Endoprosthesis Length (cm)
PLC231200*	23	12
PLC231400*	23	14
PLC271000*	27	10
PLC271200*	27	12
PLC271400*	27	14

* PLEASE NOTE: All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.

Aortic Extender Endoprosthesis

Catalogue Number	Endoprosthesis Diameter (mm)	Endoprosthesis Length (cm)
PLA230300	23	3.3
PLA260300	26	3.3
PLA280300	28.5	3.3
PLA320400	32	4.5
PLA360400	36	4.5

Iliac Extender Endoprosthesis

Catalogue Number	Distal Diameter (mm)	Endoprosthesis Length (cm)
PLL161007	10	7
PLL161207	12	7
PLL161407	14.5	7

GORE® EXCLUDER® Iliac Branch Endoprosthesis Iliac Branch Component

Catalogue Number	Diameter Iliac End (mm)	Iliac Branch Component Length (cm)
CEB231010A	10	10

Catalogue Number	Diameter Iliac End (mm)	Iliac Branch Component Length (cm)
CEB231210A	12	10
CEB231410A	14.5	10

Internal Iliac Component

Catalogue Number	Diameter Iliac End (mm)	Iliac Branch Component Length (cm)
HGB161007A	10	7
HGB161207A	12	7
HGB161407A	14.5	7

GORE® DrySeal Flex Introducer Sheath

Catalogue Number	Sheath Size (Fr)	Minimum Sheath ID (mm)	Nominal Sheath OD (mm)	Working Length (cm)
DSF1233	12	4.0	4.7	33
DSF1245	12	4.0	4.7	45
DSF1433	14	4.7	5.3	33
DSF1533	15	5.0	5.6	33
DSF1633	16	5.3	6.1	33
DSF1833	18	6.0	6.7	33
DSF2033	20	6.7	7.5	33
DSF2065	20	6.7	7.5	65
DSF2233	22	7.3	8.2	33
DSF2265	22	7.3	8.2	65
DSF2433	24	8.0	8.8	33
DSF2465	24	8.0	8.8	65
DSF2633	26	8.7	9.5	33

Catalogue Number	Sheath Size (Fr)	Minimum Sheath ID (mm)	Nominal Sheath OD (mm)	Working Length (cm)
DSF2665	26	8.7	9.5	65

GORE® Tri-Lobe Balloon Catheter

Catalogue Number	Vessel Diameter (mm)
BCM1634	16–34

GORE® Molding & Occlusion Balloon

Catalogue Number	Balloon inflation range (mm)	Balloon length (cm)	Catheter length (cm)	Wire compatibility
MOB37	10-37	4	90	0.035

INDICATIONS FOR USE: Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis Components. The GORE® EXCLUDER® AAA Endoprosthesis is intended to exclude the aneurysm from the blood circulation in patients diagnosed with infrarenal abdominal aortic aneurysm (AAA) disease and who have appropriate anatomy as described below: Adequate iliac / femoral access; Infrarenal aortic neck treatment diameter range of 19 – 32 mm and a minimum aortic neck length of 15 mm; Proximal aortic neck angulation $\leq 60^\circ$ Iliac artery treatment diameter range of 8 – 25 mm and iliac distal vessel seal zone length of at least 10 mm. Aortic Extender Endoprosthesis and Iliac Extender Endoprosthesis Components. The Aortic and Iliac Extender Endoprostheses are intended to be used after deployment of the GORE® EXCLUDER® AAA Endoprosthesis. These extensions are intended to be used when additional length and / or sealing for aneurysmal exclusion is desired.

CONTRAINDICATIONS: The GORE® EXCLUDER® AAA Endoprosthesis is contraindicated in patients with known sensitivities or allergies to the device materials and patients with a systemic infection who may be at increased risk of endovascular graft infection. **Rx Only**

**Sizing, availability, and pricing varies by country.
Please check with your Gore representative for availability.**