

Alto™

Abdominal Stent
Graft System


STENT GRAFT SYSTEM OVERVIEW

STEP 1

Select ALTO Aortic Body

Stent Graft Diameter (mm)	Aortic ID* (mm)
20	16-17
23	18-20
26	21-23
29	24-26
34	27-30


*At the intended proximal sealing ring location.



STEP 2


Select Ovation iX Iliac Limbs

Labeled Diameter (mm)	Native Iliac Vessel ID Range (mm)
10	8-9
12	10-11
14	12-13
16	14-15
18	16-17
22	18-20
28	21-25



STEP 3

Select Ovation iX Extensions



STENT GRAFT SYSTEM SPECIFICATIONS

ALTO AORTIC BODY		Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)		Stent Graft Length (mm)	
Catalog Number							
TV-AB2080-L		13	15	20		80	
TV-AB2380-L		13	15	23		80	
TV-AB2680-L		13	15	26		80	
TV-AB2980-L		13	15	29		80	
TV-AB3480-L		13	15	34		80	

OVATION IX ILIAC LIMB		Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)	Stent Graft Length with Aortic Body (mm)
Catalog Number							
TV-IL141080-I		10	12	14	10	80	125
TV-IL1410100-I		10	12	14	10	100	145
TV-IL1410120-I		10	12	14	10	120	165
TV-IL1410140-I		10	12	14	10	140	185
TV-IL1410160-I		10	12	14	10	160	205
TV-IL141280-I		10	12	14	12	80	125
TV-IL1412100-I		10	12	14	12	100	145
TV-IL1412120-I		10	12	14	12	120	165
TV-IL1412140-I		10	12	14	12	140	185
TV-IL1412160-I		10	12	14	12	160	205
TV-IL141480-I		10	12	14	14	80	125
TV-IL1414100-I		10	12	14	14	100	145
TV-IL1414120-I		10	12	14	14	120	165
TV-IL1414140-I		10	12	14	14	140	185
TV-IL1414160-I		10	12	14	14	160	205
TV-IL141680-I		11	13	14	16	80	125
TV-IL1416100-I		11	13	14	16	100	145
TV-IL1416120-I		11	13	14	16	120	165
TV-IL1416140-I		11	13	14	16	140	185
TV-IL1416160-I		11	13	14	16	160	205
TV-IL141880-I		11	13	14	18	80	125
TV-IL1418100-I		11	13	14	18	100	145
TV-IL1418120-I		11	13	14	18	120	165
TV-IL1418140-I		11	13	14	18	140	185
TV-IL1418160-I		11	13	14	18	160	205
TV-IL142280-I		12	14	14	22	80	125
TV-IL1422100-I		12	14	14	22	100	145
TV-IL1422120-I		12	14	14	22	120	165
TV-IL1422140-I		12	14	14	22	140	185
TV-IL1422160-I		12	14	14	22	160	205
TV-IL142880-I		13	15	14	28	80	125
TV-IL1428100-I		13	15	14	28	100	145
TV-IL1428120-I		13	15	14	28	120	165
TV-IL1428140-I		13	15	14	28	140	185
TV-IL1428160-I		13	15	14	28	160	205

OVATION IX ILIAC EXTENSION		Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)	
Catalog Number							
TV-EX101045-I		10	12	10	10	45	
TV-EX121245-I		10	12	12	12	45	
TV-EX141445-I		10	12	14	14	45	
TV-EX161645-I		11	13	16	16	45	
TV-EX181845-I		11	13	18	18	45	
TV-EX222245-I		12	14	22	22	45	
TV-EX282845-I		13	15	28	28	45	

CUSTOMSEAL™ POLYMER FILL KIT	Catalog Number
	TV-CS14-F

AUTOINJECTOR 2	Catalog Number
	TV-AI01-L



The ALTO™ Abdominal Stent Graft System and associated components are not available in all countries or regions.

Please contact your Endologix representative for details regarding product availability. Prior to use, refer to Instructions for Use for more information concerning Indications, Contraindications, Specific Anatomic Considerations, Warnings, Precautions, and Adverse Events.

Rx only.

©2020 Endologix Inc. All rights reserved. MM2146 Rev 03

Not for distribution in the USA.

