

Catalogue Numbers for GORE® EXCLUDER® Conformable AAA Endoprosthesis with ACTIVE CONTROL System

Trunk - Ipsilateral Leg Endoprosthesis

Catalogue Number	Intended aortic vessel diameter (mm)	Aortic Endoprosthesis diameter (mm)	Length to flow divider (cm)	Intended iliac vessel diameter (mm)	Iliac Endoprosthesis diameter (mm)	Overall device length (cm)
CXT201212	16–18	20	5.5	10–11	12	12
CXT201412	16–18	20	5.5	12–13.5	14.5	12
CXT231212	19–21	23	5.5	10–11	12	12
CXT231412	19–21	23	5.5	12–13.5	14.5	12
CXT261212	22–23	26	5.5	10–11	12	12
CXT261412	22–23	26	5.5	12–13.5	14.5	12
CXT281212	24–26	28.5	5.5	10–11	12	12
CXT281412	24–26	28.5	5.5	12–13.5	14.5	12
CXT321414	27–29	32	6.5	12–13.5	14.5	14
CXT361414	30–32	36	6.5	12–13.5	14.5	14

Aortic Extender Endoprosthesis

Catalogue Number	Intended aortic vessel diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended introducer sheath (Fr)
CXA200005	16–18	20	4.5	15
CXA230005	19–21	23	4.5	15
CXA260005	22–23	26	4.5	15
CXA280005	24–26	28.5	4.5	16
CXA320005	27–29	32	4.5	18

Catalogue Number	Intended aortic vessel diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended introducer sheath (Fr)
CXA360005	30–32	36	4.5	18

Contralateral Leg Endoprosthesis

Catalogue Number	Intended iliac vessel diameter (mm)	Iliac Endoprosthesis diameter (mm)	Contralateral Leg length (cm)	Recommended introducer sheath (Fr)
PLC121000	10–11	12	10	12
PLC121200	10–11	12	12	12
PLC121400	10–11	12	14	12
PLC141000	12–13.5	14.5	10	12
PLC141200	12–13.5	14.5	12	12
PLC141400	12–13.5	14.5	14	12
PLC161000*	13.5–14.5	16	9.5	12
PLC161200*	13.5–14.5	16	11.5	12
PLC161400*	13.5–14.5	16	13.5	12
PLC181000*	14.5–16.5	18	9.5	12
PLC181200*	14.5–16.5	18	11.5	12
PLC181400*	14.5–16.5	18	13.5	12
PLC201000*	16.5–18.5	20	9.5	12
PLC201200*	16.5–18.5	20	11.5	12
PLC201400*	16.5–18.5	20	13.5	12
PLC231000*	18.5–21.5	23	10	14
PLC231200*	18.5–21.5	23	12	14
PLC231400*	18.5–21.5	23	14	14
PLC271000*	21.5–25	27	10	15
PLC271200*	21.5–25	27	12	15
PLC271400*	21.5–25	27	14	15

* PLEASE NOTE: 16, 18, 20, 23, 27 mm Contralateral Leg Endoprostheses can be used as Iliac Extenders.

Iliac Extender Endoprosthesis

Catalogue Number	Intended iliac vessel diameter (mm)	Iliac Endoprosthesis diameter (mm)	Contralateral Leg length (cm)	Recommended introducer sheath (Fr)
PLL161007	8–9	10	7	12
PLL161207	10–11	12	7	12
PLL161407	12–13.5	14.5	7	12