

Table 24. Aortic Bifurcate Prosthesis Dimensions Sizing Guide

| Product Code | Bifurcate Main Diameter (mm) | Aortic Vessel Diameter Range Treated (mm) | Delivery System ID (F) | Delivery System OD | | Ipsilateral Length (cm) | Contralateral Length (cm) |
|--------------|------------------------------|---|------------------------|--------------------|-----|-------------------------|---------------------------|
| | | | | F | mm | | |
| AB2298 | 22 | 17.0 - 19.9 | 13 | 14 | 4.7 | 9.4 | 8.6 |
| AB2698 | 26 | 20.0 - 22.9 | 13 | 14 | 4.7 | 9.4 | 8.6 |
| AB3098 | 30 | 23.0 - 26.9 | 13 | 14 | 4.7 | 9.4 | 8.6 |
| AB3498 | 34 | 27.0 - 31.0 | 15 | 16 | 5.3 | 9.4 | 8.6 |

Table 25. Iliac Limb and Limb Extension Prosthesis Dimensions Sizing Guide

| Product Code | Limb Diameter (mm) | Iliac Vessel Diameter Range Treated (mm) | Limb Length (cm) | Delivery System OD | | Ipsilateral Treatment Length (mm) | Contralateral Treatment Length (mm) |
|--------------|--------------------|--|------------------|--------------------|-----|-----------------------------------|-------------------------------------|
| | | | | F | mm | | |
| IL1008 | 10 | 7.0 - 8.9 | 8.2 | 12 | 4.0 | 128-156 | 128-147 |
| IL1010 | 10 | 7.0 - 8.9 | 10.1 | 12 | 4.0 | 147-175 | 147-166 |
| IL1012 | 10 | 7.0 - 8.9 | 12.0 | 12 | 4.0 | 166-194 | 166-185 |
| IL1014 | 10 | 7.0 - 8.9 | 13.8 | 12 | 4.0 | 184-212 | 184-203 |
| IL1308 | 13 | 9.0 - 10.9 | 8.2 | 12 | 4.0 | 128-156 | 128-147 |
| IL1310 | 13 | 9.0 - 10.9 | 10.1 | 12 | 4.0 | 147-175 | 147-166 |
| IL1312 | 13 | 9.0 - 10.9 | 12.0 | 12 | 4.0 | 166-194 | 166-185 |
| IL1314 | 13 | 9.0 - 10.9 | 13.8 | 12 | 4.0 | 184-212 | 184-203 |
| IL1608 | 16 | 11.0 - 13.9 | 8.2 | 12 | 4.0 | 128-156 | 128-147 |
| IL1610 | 16 | 11.0 - 13.9 | 10.1 | 12 | 4.0 | 147-175 | 147-166 |
| IL1612 | 16 | 11.0 - 13.9 | 12.0 | 12 | 4.0 | 166-194 | 166-185 |
| IL1614 | 16 | 11.0 - 13.9 | 13.8 | 12 | 4.0 | 184-212 | 184-203 |
| IL2008 | 20 | 14.0 - 17.9 | 8.2 | 12 | 4.0 | 128-156 | 128-147 |
| IL2010 | 20 | 14.0 - 17.9 | 10.1 | 12 | 4.0 | 147-175 | 147-166 |
| IL2012 | 20 | 14.0 - 17.9 | 12.0 | 12 | 4.0 | 166-194 | 166-185 |
| IL2014 | 20 | 14.0 - 17.9 | 13.8 | 12 | 4.0 | 184-212 | 184-203 |
| IL2410 | 24 | 18.0 - 22.0 | 10.1 | 13 | 4.3 | 147-175 | 147-166 |
| IL2412 | 24 | 18.0 - 22.0 | 12.0 | 13 | 4.3 | 166-194 | 166-185 |
| IL2414 | 24 | 18.0 - 22.0 | 13.8 | 13 | 4.3 | 184-212 | 184-203 |

Warning: Failure to comply with sizing requirements outlined in the tables above may result in prostheses leaks, compromised flow, prostheses migration, compromised long term device durability or other complications or adverse events.

Caution: Take care to avoid excessive oversizing in tapered aortic necks.